

S-SERIES™ LEADS CASE STUDY #1

Spinal Cord Stimulation for FBSS With an S-Series Perc-Paddle Lead



D. Logé, MD
Anesthesiology and Pain Therapy Department,
General Hospital Sint-Lucas, Ghent, Belgium

Introduction and Patient History

In July 2008, a female patient presented with lumbalgia after two low back surgeries and complained of radiation pain in the groin, left leg, and foot. The patient was referred to the department of pain therapy by an orthopedic surgeon, and no paresis, paresthesia, or allodynia was found. The patient was subsequently diagnosed with failed back surgery syndrome (FBSS) with lumboschialgia, describing her pain as deep and stinging and rating her pain at 8/10 on the Visual Analogue Scale (VAS).

| | |
|-----------------------------|--------|
| Gender | Female |
| Age at time of implantation | 58 |
| Height | 1,68 m |
| Weight | 89 kg |

Table 1. Demographic data

Prior to the symptoms recurring, the patient had two weeks of pain improvement from infiltrations of the S1 nerve root. This procedure was repeated twice before the decision was made to implant her with a spinal cord stimulation (SCS) trial system. Other previous treatments included epidural injections, facet joint injections, and medication. A psychiatric evaluation revealed no mood anomalies, psychotic elements, or toxicomanic behavior, and the EEG showed a normal trace. The clinical psychological evaluation (McGill Parr Pain Questionnaire [MPQ], Pain Disability Index [PDI], Pain Catastrophizing Scale [PCS], Oswestry Low Back Pain Questionnaire, four-dimensional complaints questionnaire, Tampa scale for kinesiophobia, BECK depression scale, Utrecht Coping List [UCL], Minnesota Multiphasic Personality Inventory [MMPI]) was obtained and showed the patient to be a suitable candidate for spinal cord stimulation.

In October 2008, the patient was implanted with a 30-cm S-Series S-8 perc-paddle lead (St. Jude Medical Neuromodulation Division, Plano, TX) via percutaneous approach. A skin incision was made with an 11-blade scalpel at the L2-L3 level, after which a Touhy needle was inserted using a paramedian approach and a flat angle (<math><30^\circ</math>) to enter the epidural canal at T12/L1. The Epiducer™ lead delivery system (St. Jude Medical Neuromodulation Division, Plano, TX) was used to allow for minimally invasive placement of the perc-paddle lead. The lead was passed into the epidural canal and steered cephalad to the desired target level, with the tip of the lead at T8 and slightly to the left of the physiologic midline (Figure 1).



Figure 1. Final position of the lead with tip at T8

Intraoperative stimulation, with the Rapid Programmer™ system (St. Jude Medical Neuromodulation Division, Plano, TX), was used to confirm adequate paresthesia overlap in the patient's painful area. The patient reported paresthesia throughout the entire painful area, covering the left leg, foot, and groin as well as the left lower back area. The patient was then discharged and trialed for 28 days as required by Belgian reimbursement guidelines.

Outcome research has demonstrated that patients with FBSS benefit significantly from spinal cord stimulation.¹ The use of S-Series perc-paddle leads for percutaneous implantation and advancement in the epidural space has been found to be efficacious and safe.² In this case report, FBSS with lumboschialgia is treated with the perc-paddle lead.

S-SERIES LEADS CASE STUDY #1

Results

At follow-up after the trial (Table 2), average pain intensity decreased to 1/10 on the Visual Analogue Scale (VAS). The stimulation also resulted in a positive effect on the patient's functionality and quality of life (QoL), as measured by her rating of a number of functional activities (Tables 4 and 5) and her rating of QoL on a scale from 0–10. She also rated the stimulation as having a positive effect on her ability to perform daily activities, social activities, dependency on others, ability to perform hobbies and sports, and rest required as a result of pain (Tables 6 and 7). The stimulation was rated excellent in its ability to help with her pain problem, and she also indicated that she would definitely redo the procedure (Tables 8 and 9). In addition, compared to baseline, the patient's quality of sleep and problems falling asleep and waking up during the night due to pain improved during the trial period (Table 7).

After the successful trial, the lead was connected via a 60-cm extension to a GenesisXP™ primary cell implantable pulse generator (St. Jude Medical Neuromodulation Division, Plano, TX) and implanted in a subcutaneous abdominal pocket. Parameters programmed were as follows: electrode polarity: 4+ 5- 6-; frequency: 30 Hz; pulse width: 450 µs; perception amplitude: 2,5 mA; comfort amplitude: 4,0 mA; Maximum Tolerable amplitude: 6,5 mA.

| | VAS Before | VAS After |
|------------------------------|------------|-----------|
| Current pain (at visit) | 2,9 | 0,0 |
| Worst pain in the last week | 10,0 | 5,9 |
| Lowest pain in the last week | 1,9 | 0,0 |
| Mean pain in the last week | 5,1 | 1,0 |

Table 2. Pre- and postoperative pain scores

| | Before | After |
|-------------------|--------|-------|
| Extremely intense | 30% | |
| Intense | | |
| Strong | 30% | |
| Mild | 40% | 5% |
| Light | | |
| No pain | | 95% |

Table 3. Description of the pain variation a day before and after implantation

| | Without problems | Small problems | Some problems | A lot of problems | Impossible |
|--------------------------------|------------------|----------------|---------------|-------------------|------------|
| Lift heavy weights (+10 kg) | | | | | X |
| Sitting down 30 min | | X | | | |
| Walk 500 m | | | | | X |
| In and out of bed without help | | X | | | |
| Dress without help | | X | | | |
| Take 10 steps on stairs | | | | | X |
| Shower/bathe without help | | | X | | |

Table 4. Activity questionnaire according to Belgian Reimbursement requirements at baseline

| | Without problems | Small problems | Some problems | A lot of problems | Impossible |
|--------------------------------|------------------|----------------|---------------|-------------------|------------|
| Lift heavy weights (+10 kg) | | | X | | |
| Sitting down 30 min | X | | | | |
| Walk 500 m | X | | | | |
| In and out of bed without help | X | | | | |
| Dress without help | X | | | | |
| Take 10 steps on stairs | X | | | | |
| Shower/bathe without help | X | | | | |

Table 5. Activity questionnaire according to Belgian Reimbursement requirements after implantation

| | Before | After |
|--|--------|-------|
| Daily activities (cleaning, cooking, shopping, etc.) | 6,3 | 1,3 |
| Social activities (going out, meet with friends, etc.) | 4,6 | 1,7 |
| Dependency of other people | 4,6 | 0,1 |
| Hobbies, sports, leisure | 9,3 | 5,7 |
| Rest due to the pain (bed/couch) | 3,6 | 0,1 |

Table 6. Influence of the pain on daily activities and rest in the week, measured before and after implantation. 0 = Pain does not inhibit the activity at all; 10 = Pain totally prevents doing the activity. 0 = No rest necessary due to the pain; 10 = Rest the whole day.

| | Before | After |
|---------------------------|---|---|
| Quality of sleep | Very good Good Mild Bad Very bad | Very good Good Mild Bad Very bad |
| Problems falling asleep | Never Some nights Most nights Every night | Never Some nights Most nights Every night |
| Waking up due to the pain | Never Some nights Most nights Every night | Never Some nights Most nights Every night |

Table 7. Answers to questions on sleeping patterns

Did the therapy help with your pain problem?

| X | Excellent | | Little |
|---|-----------|--|-----------------------|
| | Very good | | Not at all |
| | Good | | Increased the problem |
| | Mild | | |

Table 8. Answer to the question “Did the therapy help with your pain problem?”

Would you redo the procedure?

| X | Definitely yes |
|---|----------------|
| | Yes |
| | Maybe |
| | No |
| | Definitely not |

Table 9. Answer to the question “Would you redo the procedure?”

Conclusion



Figure 2. Lateral view on the lead to confirm posterior epidural lead placement. The trial cable connection can also be seen on the lower left of the picture.

This case report describes a successful outcome using a novel approach to implant S-series perc-paddle leads in a patient diagnosed with FBSS and lumboschialgia. Following implant, the patient had three control visits at which no reprogramming was necessary and paresthesia was reported in the entire painful area, including the low back. In August 2009 (10-month follow-up), the patient’s VAS score remained 1/10. Tramadol was entirely stopped and paracetamol was only used as needed. Our result, as reported in this case report, has been encouraging.

One key to this successful outcome was performing intraoperative testing of the lead placement while the patient was awake (Figure 2). The case suggests that a perc-paddle lead implanted with a novel percutaneous implant technique is safe and efficacious for failed back surgery syndrome, as has been reported previously.³

All data, including VAS, pain characteristics, medication intake, functional activities, and QoL parameters were collected using the Belgian Pain Society evaluation form. This form was developed by the Belgian Pain Society in 1997 and is ruled necessary by the Belgian government for reimbursement of all patients with chronic pain treated with implantable devices. According to Belgian law, an application for Institutional Review Board (IRB) was not required (at the initiation of this case report) for studies that concern regular medical practice and follow-up, and thus no application for IRB approval was made.

Our result, as reported in this case report, has been encouraging. One key to this successful outcome was performing intraoperative testing of the lead placement while the patient was awake. The case suggests that a perc-paddle lead implanted with a novel percutaneous implant technique is safe and efficacious for failed back surgery syndrome as has been reported previously.

— D. Loge, MD³

St. Jude Medical is focused on reducing risk by continuously finding ways to put more control into the hands of those who save and enhance lives.

1. Mekhail NA et al. Clinical applications of neurostimulation: forty years later. *Pain Practice*. 2010; 10(2): 103-122.
2. Logé D, De Coster O, Washburn S. A retrospective data collection study to evaluate the feasibility and safety of percutaneous introduction of a narrow paddle lead into the epidural space. Poster presented at: WIP Congress; March 2009; New York, NY.
3. Logé D, De Coster O. A prospective evaluation to assess the safety of using a newly developed delivery system for percutaneous introduction of SCS paddle leads into the epidural space. Poster presented at: Annual meeting of the North American Neuromodulation Society; Dec 3-6, 2009; Las Vegas, NV.

ATRIAL FIBRILLATION

CARDIAC RHYTHM MANAGEMENT

CARDIOVASCULAR

NEUROMODULATION

Global Headquarters

One St. Jude Medical Drive
St. Paul, Minnesota 55117
USA
+1 651 756 2000
+1 651 756 3310 Fax

Neuromodulation Division

6901 Preston Road
Plano, Texas 75024
USA
+1 972 309 8000
+1 972 309 8150 Fax

St. Jude Medical Europe, Inc.

The Corporate Village
Da Vincilaan, 11 Box F1
1935 Zaventem
Belgium
+32 2 774 68 11
+32 2 772 83 84 Fax

sjmneuro.com



ST. JUDE MEDICAL™

MORE CONTROL. LESS RISK.

Not all products approved in all countries. Please ask your local representative.

St. Jude Medical SCS devices are intended to aid in the management of chronic intractable pain of the trunk and/or limbs, which may have derived from specific disease states or diagnoses. St. Jude Medical neurostimulation systems are not intended to treat or cure specific disease states or diagnoses.

Indications for Use: Spinal cord stimulation as an aid in the management of chronic, intractable pain of the trunk and limbs. **Contraindications:** Demand-type cardiac pacemakers, patients who are unable to operate the system or who fail to receive effective pain relief during trial stimulation. **Warnings/Precautions:** Diathermy therapy, cardioverter defibrillators, magnetic resonance imaging (MRI), explosive or flammable gases, theft detectors and metal screening devices, lead movement, operation of machinery and equipment, postural changes, pediatric use, pregnancy, and case damage. Patients who are poor surgical risks, with multiple illnesses, or with active general infections should not be implanted. **Adverse Events:** Painful stimulation, loss of pain relief, surgical risks (e.g., paralysis). Clinician's manual must be reviewed prior to use for detailed disclosure. Rx only.

Epiducer, GenesisXP, Rapid Programmer, and S-Series are trademarks of Advanced Neuromodulation Systems, Inc. d/b/a St. Jude Medical Neuromodulation Division. ST. JUDE MEDICAL, the nine-squares symbol, and MORE CONTROL. LESS RISK. are trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2011 St. Jude Medical Neuromodulation Division. All rights reserved.

Item 0782-01

